


# **Stainless Steel Filter Holders**

Vacuum filtration holders made from nonbreakable stainless steel are ideal for applications where durable and sterilizable equipment is required.

This stainless steel unit may be flame sterilized to sanitize the holder between samplings. The Set pins and locking nut on the funnel help prevent twisting and tearing of the membrane. This holder has the added benefit of being autoclavable with a membrane in place.

### **SPECIFICATIONS**

KSF 47	352600	
Materials	Type 304 stainless steel funnel, base, and support screen	
	Silicone rubber stopper	
	PTFE gasket	
Filter Size	47 mm	
Prefilter size	35 mm	
Filtration area	9.6 cm2	
Connections	#8 stopper, mounts in a standard 1L filtering flask or standard cup of any Sterlitech filtering manifold	
Pressure	21 in Hg vacuum pressure recommended	



#### **CLEANING AND MAINTENANCE**

- Routinely inspect surfaces for nicks and rough edges.
- Clean all components with a standard, low-alkaline, non-abrasive detergent such as Alconox or Liquinox.
- Rinse with hot tap water followed by a final rinse with distilled or deionized water.
- Store dry.
- Avoid acids.
- To autoclave, wrap assembly in Kraft paper or other steam permeable paper. Autoclave for 15 minutes at 121°C (250°F); SLOW exhaust. Holders can be used after they have cooled and can be handled without gloves.

#### **ASSEMBLY AND OPERATION**

- 1. Holding base upright and using Tweezers (Pad) Miracle Tip Ends or blunt-tipped forceps:
  - a. Place Teflon gasket on base.
  - b. Place support screen on top of gasket.
  - c. Center membrane filter on support screen.
  - d. Center prefilter (if used) on top of membrane.
- 2. Using set pins to center, place funnel on top of base assembly.
- 3. Twist locking nut to tighten.
- 4. Wet stopper with water to lubricate then slide stopper onto base.
- 5. Mount in flask or manifold.

## **ORDERING GUIDE**

Funnel	352601
SS support screen	311504
PTFE gasket	311505
Base	312602
#8 stopper	311404