UPGRADE of Stainless Steel Filter Holder

SF, Stainless Steel Filter Holder



SS, Stainless Steel Filter Holder



SF 1	SF 3	SF 5	Model	SS 1	SS 3	SS 5
180100-10	180100-30	180100-50	Cat. No.	167103-10	167103-30	167103-50
100 ml	300 ml	500 ml	Capacity	100 ml	300 ml	500 ml
SS316			Material	SS316		
Outer Spin-Lock			Connection	Inner Spin-Lock		
Yes			Membrane Guide 1	No		
Yes			Forceps Access Points 2	No		
Laser-Marked			Funnel Graduation 3	Engraved-Ring		
11.6 cm ²			Effective Filtration Area	9.6 cm ²		
• Rocker, SF Funnel			Compatible with	Rocker, SS Funnel		
• Pall, 4870/4871 Disposable Funnel						



- Membrane guide helps position membrane in the middle of support base easily.
- 2 Forceps access points help retrieve membrane easily without any damage.
- 3 Laser-Marked graduation helps add sample accurately with less residue.



